


Unison UK-C8

Polar to Non-Polar Compounds, Basic Compounds

Better-balancer and higher-speed separation than ODS



Conventional C8 columns are known for short retention times and faster elution than their ODS counterparts. Unison UK-C8 provides optimal hydrophobicity and electrostatic interaction - enabling scientists to reduce analysis time and achieve acceptable separation.

Analytical Columns

Semi-Prep Columns

Unison UK-C8 3 µm Particle Size Stationary Phase	Length mm	Internal Diameter				
		1	2	3	4.6	6
	10		UK820	UK830	UK800	
20		UK829	UK839	UK809		
30	UK811	UK821	UK831	UK801	UK861	UK8P1
50	UK812	UK822	UK832	UK802	UK862	UK8P2
75	UK813	UK823	UK833	UK803	UK863	UK8P3
100	UK814	UK824	UK834	UK804	UK864	UK8P4
150	UK815	UK825	UK835	UK805	UK865	UK8P5
250	UK816	UK826	UK836	UK806	UK866	UK8P6
500				UK807		

Guard Cartridges

Size	Code	Pieces
1mm	GCUK8C	3
2-6mm	GCUK8S	3
10mm	GCUK8M	2

Guard Holders

Column Coupler Included	
1-6mm	GCH01S
10mm	GCH02M